

# NexScope

NEPHROSCOPE AND MINI NEPHROSCOPE



COMPETENCE IN SURGICAL LASER SINCE 1989  
LASER TECHNOLOGY - MADE IN GERMANY  
[WWW.LISALASER.DE](http://WWW.LISALASER.DE)



## KEY FEATURES

- ★ Outstanding performance: contrast enhanced image bundle, best vision
- ★ Safe and reliable: robust, tightly sealed, laser welded, easy to use and to clean
- ★ Direction of view ideally aligned at 8°
- ★ Strong irrigation flow enables high stone-free rates
- ★ Long lifetime for lowest cost of ownership
- ★ Lightweight design
- ★ Detachable instrument bridge for optimized cleaning and inspection
- ★ Compatible with most popular laser and lithotripter systems
- ★ Autoclavable: 5 min. @ 134 °C
- ★ Made in Germany



## Technical Specifications **NexScope**

### Nephroscope and Mini Nephrosopes

<b>NexScope</b>	<b>Description</b>	<b>Part number</b>
NexScope 19 F 220	working channel 10 F, field of view 84°	101 503 437
NexScope Sheath 24 F	with obturator, colour code: red	101 503 440
NexScope Cleaning Brush, Working Channel 10 F	colour code: blue	101 503 486

### **Mini NexScope**

Mini NexScope 13.5 F 220	working channel 6 F, field of view 84°	101 503 435
Mini NexScope 13.5 F 320	working channel 6 F, field of view 84°	101 503 436
NexScope Sheath 18 F	with obturator, colour code: yellow	101 503 439
NexScope Cleaning Brush, Working Channel 6 F	colour code: white	101 503 487

### **Consumables and accessories for NexScope and Mini NexScope**

NexScope Instrument bridge	for <b>NexScope</b> and <b>Mini NexScope</b>	101 503 438
Cleaning Brush for Side Port	of instrument bridge, colour code: white	101 503 488
Cleaning Brush for Main Port	of instrument bridge, colour code: blue	101 503 489
Duck Bill Valves	disposable, 1 package á 10 pieces	101 503 498
Sealing caps, for fibre fixation	disposable, 1 package á 10 pieces	101 503 499
Stop cock plugs, single use	1 package á 10 pieces	101 503 509

CE acc. Council Directive 93/42/EEC

U.S. federal law restricts these devices to sale by or on the order of a physician

Specifications are subject to change without notice  
Made in Germany 2017 - 11  
Brochure Part No. 036 013 014



LISA laser products OHG  
Albert-Einstein-Str. 1 - 9  
37191 Katlenburg-Lindau  
Germany  
fon: +49 5556 9938-0  
fax: +49 5556 9938-10  
info@lisalaser.de www.lisalaser.de



**LISA LASER PRODUCTS**  
COMPETENCE IN SURGICAL LASER

